

Endometrial Cancer

- 1. SCOUT**
- 2. CORONAL HASTE** – FOV to include kidneys and bladder
- 3. CORONAL T2 HASTE FS THICK SL** – FOV to include kidneys and bladder
- 4. AXIAL T2** – entire pelvis
- 5. SAGITTAL T2 FSE** – small FOV, high resolution
- 6. SAGITTAL T2 BLADE** – small FOV, high resolution
- Decide FSE vs. BLADE for remainder of T2WI, MD check if need help
- 7. CORONAL OBLIQUE T2** – along the long & short axis of the uterus
- 8. SAGITTAL DIFFUSION** – b values of 50, 400 and 800
- 9. CORONAL OBLIQUE (short axis) DIFFUSION** – b values of 50, 400 and 800
- 10. AXIAL T1**
- 11. AXIAL VIBE** – entire pelvis

MD check for plane of VIBEs

- 12. SAGITTAL and CORONAL OBLIQUE (short axis) VIBE**

Inject contrast

- 13. DYNAMIC VIBEs** – same plane as pre-VIBE starting with sagittal and alternating with short axis, to 4 minutes, with subtractions
- 14. AXIAL VIBE** – entire pelvis