

Liver Transplant

1. SCOUT
2. CORONAL TRUEFISP BH LOCALIZER
3. AXIAL T2 HASTE MBH
4. CORONAL T2 HASTE MBH
5. AXIAL IN and OPPOSED PHASE T1 GRE BH
6. AXIAL DIFFUSION – b values of 50, 400, 800
7. CORONAL LOCALIZER non-BH
8. AXIAL T2 BLADE FS resp nav
9. T2 SPC MRCP non-BH – navigator
10. HASTE MRCP – radial stack
11. AXIAL VIBE

Inject contrast

12. MRA
13. AXIAL VIBE – 70 sec.
14. CORONAL VIBE – 2 min.
15. AXIAL VIBE – 3 min.
16. AXIAL VIBE – 5 min.

OPTIONAL:

T2 HASTE FS THICK SL PAT2

AXIAL AND CORONAL T2 HASTE BH – thin 3mm, cover pancreatic duct & CBD