

Liver - Eovist

- If MRCP needed, it must be performed pre-contrast
- Use for suspected FNH, do not use in cirrhotic patient or with hemangiomas

1. SCOUT

2. CORONAL TRUEFISP BH LOCALIZER

3. AXIAL IN and OPPOSED PHASE T1 GRE BH

4. AXIAL DIFFUSION – b values of 50, 400, 800

5. AXIAL VIBE BH

Inject Eovist contrast (0.0025 mmol/kg rounded up to the full bottle)

6. AXIAL VIBE – triple arterial phase

7. AXIAL VIBE – 60 sec.

8. CORONAL VIBE – 90 sec.

9. AXIAL VIBE – 3 min.

10. AXIAL T2 HASTE MBH

11. AXIAL VIBE – 5 min.

12. CORONAL T2 HASTE MBH

13. AXIAL VIBE – 10 min.

14. CORONAL LOCALIZER non-BH

15. AXIAL T2 BLADE FS resp nav

16. AXIAL VIBE – 15 min.

Access our Modality Protocols site for all current protocols at:
www.MIARAD.com

INTEGRITY | TEAMWORK | EXCELLENCE | SERVICE | ACCOUNTABILITY

Reviewed 6-27-24